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How to fabricate a CA<sup>®</sup> CLEAR ALIGNER lower splint  
using the VECTOR<sup>®</sup> 40 expansion screw acc. to Prof. Hinz

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## How to fabricate a CA® CLEAR ALIGNER lower splint using a VECTOR® 40 expansion screw acc. to Prof. Hinz

Special miniaturised screw with black directional arrow for fabrication of CA® splints, allowing expansions in the lower jaw (up to 3 mm). The retention arms allow safe anchorage using cold curing resin.

By incorporating specific VECTOR® expansion screws in CA® CLEAR ALIGNER splints, treatment of mal-occlusions such as crowdings or rotations of single teeth or groups of teeth gets even quicker and more efficient and economical in comparison with the conventional CA® set-up technique. Aesthetic qualities and transparency of the splints remain unchanged.

## We recommend the following specialized literature for further information:

CA® Book Volume 2 Advanced  
English: REF 5309 / German: REF 5307

## Important:

Make sure not to exceed the maximal tooth movement in the CA® CLEAR ALIGNER technique acc. to Dr. Pablo Echarri (CA® Book volume 2, Advanced, pages 52/53).

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Skalpell

**Working material for pressure moulding:**

REF 2487 VECTOR® 40  
REF 3207 ISOFOLAN®  
REF 3407 CA® Foil hard (0.75mm)  
REF 3218 BIOCRYL®-Resin monomer  
REF 3221 BIOCRYL®-Resin polymer  
REF 3701 SD Cyano Veneer Fast  
REF 3703 SD Pipette

Pressure moulding device

**Additionally recommended material:**

REF 3442 SIL-KITT red  
REF 5288 Diamond separating disc  
REF 2041 CARAT® Flat pliers  
REF 5372 Set-up manual saw  
REF 3606 SD Separating disc, serrated  
REF 3265 LC bur C  
REF 3460 SD Foil scissors A  
REF 3461 SD Foil scissors B  
REF 3226 BIOPERM® Trimmer  
REF 3382 DIMO® brown, coarse  
REF 3381 DIMO®PRO  
REF 3710 SD Polishing paste A  
REF 3711 SD Polishing paste B  
REF 3712 SD Polishing brush  
REF 3713 SD Felt wheel  
REF 3714 SD High gloss polishing disc  
REF 3372 POLYFIX disc

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Block out the lower jaw model in the area of crowdings and undercuts using SIL-KITT. Mark the midline.

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
Place the model on the model platform. Thermoform ISOFOLAN® according to the manufacturer's instructions, reduce the model at the model rim and make interdental incisions.

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Remove VECTOR® 40 wax space holder at the surface.

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Fix VECTOR® 40 expansion screw with a drop of CA® CYANO VENEER FAST on ISOFOLAN® prior to positioning.

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If necessary, adapt guiding pins to jaw width.

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**Note:**

Reducing the guiding pins affects expansion width!

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The dental arcade is placed on the model platform preferably with the front facing the locking shaft. (Full models should be embedded in the granules).

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Fix a CA® foil hard and mix BIOCRYL® Resin until it has a creamy consistency. Apply a few drops of BIOCRYL® Resin to the metal parts during the heating time of the CA® foil. During the last 5 seconds of the heating time more resin should be added generously. Once the heating time has ended, close the pressure chamber and thermoform.

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Ignore any acoustic signals during this period, making sure not to evacuate or to switch off the machine.

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
**Important:**

The model should stay under pressure in the unit for at least 10 minutes.

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Roughly reduce material using SD Foil scissors B.

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Cut material along the midline, with the splint remaining on the model for stability reasons.

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
The area around the screw is carefully thinned out using a bur.

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Expose the screw spindle using a scalpel.

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Reduce material to required length and separate it carefully from the lingual side.

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
Expose the transition to the screw using a scalpel.

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
Remove wax residues from the bore holes using an activation key.

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Finish and polish the splint as usual.

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We recommend using the material listed on page 3.

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

**Hint:**

Please note this instruction does not replace attending our technical CA® courses.

We also recommend participating in the following additional course:

„Course No. 17 - module 5: CA® CLEAR ALIGNER splints with specific VECTOR® expansion screws“.

For further information and details on the course please visit our website:

→ [www.scheu-dental.com](http://www.scheu-dental.com) → SCHEU-ACADEMY → Seminar → Topics: CA® CLEAR ALIGNER

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More instructions on how to fabricate aligners and splints with different VECTOR® screws:

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REF 2484, REF 5317

**CA® CLEAR ALIGNER  
with VECTOR® 30 screws  
and CA® Arch CHROMIUM**

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REF 2487

**CA® CLEAR ALIGNER  
with VECTOR® 45 screw  
acc. to Prof. Hinz**

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REF 5329

**TFA Splint  
with TFA Kit  
acc. to Dr. Echarri**

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REF 5319

**CA® CLEAR ALIGNER  
with CA® Distalizer Kit  
acc. to Dr. Echarri**