IMPRIMO® LC Splint flex

Produktinformation

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IMPRIMO® LC Splint flex

Product information

Description:	IMPRIMO® LC Splint flex is a light-curing methacrylate-based resin that is specifically engineered for production of high-precision transparent occlusal splints as well as for bruxism and lower jaw protrusion splints. Photo-polymerization is done by means of image projection systems (385 nm). IMPRIMO® LC Splint flex has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Splint flex should be used in building processes with a z-resolution of 50 - 100 microns. A minimum wall thickness of 1.5 mm is required for reasons of stability. Printer parameters available for Asiga MAX™ and Asiga PRO 4K.
Processing instructions:	Make sure IMPRIMO® LC Splint flex is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.
Contra-indications:	IMPRIMO® LC Splint flex should not be used for any purpose other than the generative CAD / CAM fabrication of high-precision, transparent occlusal splints as well as of bruxism and lower jaw protrusion splints. IMPRIMO® LC Splint flex is not suitable for printing aligners. IMPRIMO® LC Splint flex is a medical device class I.
Post-processing:	Upon completion of the printing process, the objects shall be removed from the model platform. Uncured resin residues can be removed completely using Isopropanol and a soft brush or can be resolved in the cleaning device IMPRIMO® Clean using IMPRIMO® Cleaning Liquid. Dry parts, until there are no residues of isopropyl alcohol (IPA) or equivalent cleaning agent.
Post exposure:	Complete curing of the objects is done in the curing unit IMPRIMO® Cure. The use of protective gas to prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. The polymerization process in the light oven takes 2 x 10 minutes in a nitrogen environment. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.
Shelf life / storage:	Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.
Hazard warnings:	May cause allergic skin reaction. May be harmful to aquatic organisms with long lasting effects. Attention: Polymerized resins are chemically resistant, avoid stains on clothes.
Disposal:	According to local authority prescriptions. Cured objects can be disposed of with other household waste. Un-polymerized material has to be disposed of at a waste collection point. All instructions on the processing of our materials - in oral, written or practical form - are given to the best of our knowledge and must be read as references. Material use and processing is beyond our control and lies within the responsibility of the user. For professional use only.
Technical data:	Density: 1,1 g/L Viscosity: < 1,5 Pa s Shore hardness D: ≥ 65 Flexural strength: ≥ 2 MPa Flexural modulus: ≥ 100 MPa Water absorption: ≤ 32 µg/mm ³ Impact resistance: ≥ 50 J/m Elongation at break: ≥ 40 % Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity) *These data result from measurements of a representative sample. Contains: Alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers, dyes and pigments. Material Safety Data Sheet: www.scheu-dental.com/en/service