

Quad Helix



Important information!

The Quad Helix should be securely ligated into the lingual/palatal sheath with steel or elastic ligatures!

FOR SINGLE USE ONLY

SCHEU-DENTAL and his agents will not be responsible for any or all claims resulting from multiple use of this product.

WARNING

This product contains nickel and chromium. Patients known to be allergic to these metals should not be treated with this product.

CE 0044

SCHEU-DENTAL GmbH
www.scheu-dental.com

Am Burgberg 20
58642 Iserlohn · Germany

SCHEU
Dental Technology

DE/GB 500/11/19 G REF PM0115.01

