

DEXIS™ Intraoral Scanner Family System Requirements

The IS 3800W, IS 3800 and IS 3700 are intraoral scanning devices that are used in combination with IS ScanFlow scanning acquisition software to acquire digital impressions. It is essential that the computer meets the system requirements to ensure optimal performance of your intraoral scanner.

Computer System Requirements

If necessary, you must adapt your computer system configuration to meet the requirements below:



IMPORTANT

It is MANDATORY to check that your system configuration is compatible with the Computer System requirements for IS ScanFlow software.



NOTE

Always use Microsoft® Windows® Update to ensure the latest security patches are correctly installed.

The products, equipment and services illustrated and described in this brochure reflect knowledge at the time of printing. Dental Imaging Technologies Corporation accepts no liability for any deviation from the illustrations in terms of color or shape, or any errors or print errors, and retains the right to make changes to the brochure at any time. Full or partial reprinting is only permitted with permission from Dental Imaging Technologies Corporation. For indications for use, please visit: DEXIS.com

Dental Imaging Technologies Corporation reserves the right to make changes to specifications and features shown herein, or to discontinue the product described at any time without notice or obligation. Contact your local authorized representative for the most current information. Manufactured by Dental Imaging Technologies Corporation, 450 Commerce Drive, Quakertown, PA 18951.

© 2024 Dental Imaging Technologies Corporation.
All rights reserved. DX00118\Rev B
All trademarks are the property of their respective owners.

Item	Recommended
CPU	Intel Core i7, 9th generation
RAM	32 GB RAM
Monitor	Screen resolution of 1920 X 1080
Operating system	Windows 10 Professional version 1809 or higher
USB port	IS 3800W: USB 3.0 for WiFi adapter IS 3800: USB Type-C IS 3700: USB Type-A 2.0 high speed
Video card	NVIDIA GeForce RTX 2060, 6 GB memory
	NVIDIA GeForce RTX 3060, 6 GB memory for Al Matching
Video card driver	Support OpenGL 4.3 and OpenCL 1.1

The computer and its screen should be situated in or close to the operating areas where the practitioner can easily view it when using the intraoral scanner.



DEXIS | 450 Commerce Drive Quakertown PA 18951 | USA 1-888-88-DEXIS