



YOUR RETAINER WORKFLOW JUST TURNED CHAIRSIDE

Did you know that you can produce same day in-house retainers in just 20 minutes?

Learn all about it at our Digital Workflow event.

Thursday 7th November
6:30pm - 8:30pm
MAGNA Science Adventure Centre
S60 1DX Rotherham



A 3D printer can be the starting point for a frictionless chairside retainer workflow.

3D printing unlocks the full potential of your MINISTAR S® and 3DISC Heron intraoral scanner for your orthodontic treatment. Paired with the dedicated SprintRay Chairside 3D Printing portfolio, all you need to do is print your scan, thermoform your retainer and make your patients happy.

Register today at
sales@eurodentic.co.uk



3DISC HERON IOS DESIGN SERVICES PRO 95 3D PRINTER PRO 55 3D PRINTER PROCURE 2 PROWASH S SPRINTRAY MATERIALS MINISTAR S

E-MAIL: SALES@EURODONTIC.CO.UK
TEL: +44 (0) 114 2766813